



# Bien 155: Biotechnology Lab

## Thursday Lab Section 021

Group 1

pET28(b)-CyPet-SUMO1

Amber Alden, Dev Bhatt, Chulin Cao, Viraj Khatri



# Introduction

## CyPet

- Cyan fluorescent protein
- derived from *Aequorea victoria* (also called: crystal jelly)

## SUMO1 Gene

- Ubiquitin-like modifier, small
- SUMOylation

## Top10 Electrocomp™ E.coli (Invitrogen) bacterial cells


- E.coli is not naturally competent-> electroporation
- Fast: No  $lacI^q$  gene

## BL<sub>21</sub> cells


- Transforms plasmid DNA into protein expressed bacterial cells
- Isopropyl- $\beta$ -D-thio-galactoside (IPTG)
- T7 RNA polymerase



**Figure 1:** SUMO1 gene ribbon structure



**Figure 2:** Lab 1 cell culture results




**Figure 3:** Lab 4 transformation plate results

# Introduction

## Bradford Assay (Coomassie Blue Dye)


- Colorimetric method of protein quantitation (concentration)
- Spectral shift
  - Free dye:
    - Abs max at 470 nm (brown color)
  - Protein bound:
    - Abs max at 595 nm (blue color)
- Coomassie dye ligands
  - Proportional amount to positive charges on protein
  - Bind to basic amino acids (arginine, lysine)



<sup>2</sup>Figure 4: Bradford Assay

## Förster Resonance Energy Transfer (FRET) Assay

- [Donor to acceptor] Non-radiative energy transfer phenomenon
- Emission and excitation spectrum-overlap significantly
  - Fluorescent emission always has a longer wavelength (lower energy state) than excitation



Clegg, R.M., *Curr Opin Biotechnol* (1995)

Sapsford, K.E *Angew Chem Int Ed Engl* (2006)

Figure 5: Overview of FRET Assay





Figure 6: FRET Assay Analysis



# Experimental Design Timeline



# Experimental Procedure

## Week 1 -Week 3 CyPet-SUMO1 plasmid DNA preparation


- Mix competent cells with plasmid DNA and perform eletroporation.
- Incubate the cells on LB agar plate overnight then pick up colonies into medium with Kanamycin.
- Lyse bacteria and clear it by centrifuge then the plasmid is ready to purification.
- QIAprep membrane only absorb DNA after washing and elution the high-quality plasmid is made.
- Check the gene clones by digestion and PCR.
- Retrieve sequences of SUMO1 from NCBI and assemble it with the sequence of CyPet.
- Use CLC Sequence Viewer and SnapGene/4peak to verify the sequence.



# Experimental Procedure

## Week 4 -Week 8 **Protein Extraction and determination**


- Perform electroporation with BL21 cells and CyPet-SUMO1 plasmid DNA mixture.
  - Incubate the mixture on LB agar plate with kanamycin.
  - Lysing cells by sonication 5 sec on and 5 sec off intervals.
  - Centrifuge to separate DNA and proteins.
  - Pour the target DNA solution into column with washed Ni-NTA beads inside.
  - Wash the protein with 3 different buffers.
  - Inject the protein into dialysis bag after elute it.
  - Put the dialysis bag in PBS overnight.
- 
- Prepare BSA standards and make standard curve on laptop.
  - Calculate the concentration of our sample and perform electrophoresis with 5 microgram of sample.
  - Check the fluorescent emission intensities of CyPet at 475 nm and use the given standard curve to identify the concentration of the fluorescent proteins.
  - Compare CyPet-SUMO1 with YPet-Ubc9 by FRET assay.



# Results

## Digestion Checking:

Cuts DNA with restriction enzymes




Cell Culture




Mutation in DNA sequence

## Digest checking 1




## PCR checking 1




PCR Checking:  
BP band in the PCR run



# Results




Week 4: LB Agar Plate after incubation



Week 6: Protein transferred into dialysis bag for being dialyzed overnight

Abs vs. Conc




BSA standard to measure concentration of unknown protein and total protein concentration was 1932 μg/mL




# Results

## Week 8: Analysis of produced proteins via FRET Assay


Lab 7 Bradford Coomassie Assay Absorbance Data



FRET Spectrum



1μM CyPet SUMO1 -10μM YPet-UBC9



$$\text{Fluorescence Protein Concentration} = \frac{63979.469}{4370.9} = 14.63\mu\text{M}$$

$$\text{Bradford Essay Concentration Conversion: } \frac{1.932\text{g/L}}{38000\text{ g/mol}} = 50.84\mu\text{M}$$

$$\text{Protein Purity Ratio} = \frac{14.63\mu\text{M}}{50.84\mu\text{M}} \times 100 = 29.1\%$$

$$\text{FRET Signal} = \frac{\text{Emission 530/414}}{\text{Emission 475/414}} = \frac{3211.78}{4574.59} = .70209$$

# Conclusion

## Labs 1-3

- Cloned CyPet-SUMO1 via E.Coli transformation and selection
- Used sequenced data to identify errors in sequence as device error or genuine mutations
- Performed Agarose Gel Electrophoresis to visualize DNA bands and verify cloning progress

## Labs 4-6

- Transformed CyPet-SUMO1 plasmid DNA into BL21(DE3) E.Coli cells for protein expression
- Extracted recombinant proteins and utilized Ni-His tag purification

## Labs 7-8

- Quantified yield using Bradford Assay and FRET
- Obtained purity ratio (29.1%) and FRET signal (.70207)


# Future Directions

## FRET:


- Potential use in measuring and building protein charge transfer systems
- Analytical method in construction of large synthetic protein structures

## Post translational modification systems:

- Potential use in Adenoviral Vector and Lentiviral vector manufacturing
  - Emergence of viral therapeutics



Source: Signa Laboratories



Belcher et al. , MIT BE

<https://doi.org/10.1126/science.1122716>

# References

A. (2021). *WHICH PROTEIN ASSAY IS BEST FOR YOU?* azure biosystems. <https://azurebiosystems.com/blog/protein-assay-best/>

<sup>1</sup>Alonso, A., Greenlee, M., & Matts, J. (2022). *Emerging roles of sumoylation in the regulation of actin, microtubules, intermediate filaments, and septins*. ResearchGate. [https://www.researchgate.net/figure/Structure-of-ubiquitin-and-SUMO-proteins-Ribbon-drawing-of-ubiquitin-Smt3-SUMO1-and\\_fig3\\_277414289](https://www.researchgate.net/figure/Structure-of-ubiquitin-and-SUMO-proteins-Ribbon-drawing-of-ubiquitin-Smt3-SUMO1-and_fig3_277414289)

(2022). *SUMO1 small ubiquitin like modifier 1 [ Homo sapiens (human) ]*. NIH: National Library of Medicine. <https://www.ncbi.nlm.nih.gov/gene/7341>

<sup>2</sup>(n.d.). *BBS OER LAB MANUAL*. OPEN LIBRARY. <https://ecampusontario.pressbooks.pub/biochem2l06/chapter/3-1-lab-overview-and-background-information/>

Liao, J (n.d). *Biotechnology Lab (Bien 155L) Manual*: University of California, Riverside

Signagen Laboratories , Introduction of Adenovirus <https://signagen.com/Introduction-to-Adenovirus>

Virus-Enabled Synthesis and Assembly of Nanowires for Lithium Ion Battery Electrodes: <https://www.science.org/doi/full/10.1126/science.1122716>